# OPTIMAX

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*Perfect microprocessor control allows easy setting of the developer temperature.*

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
*Easy to service high-precision racks guarantee a gentle film advance and makes daily care easy.*

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*The special combination of rollers in the dryer rack and the fan ensure gentle drying of the films.*

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*OPTIMAX: The compact allrounder that saves you money.*

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**Quality leads to success:  
worldwide.**

#### **Multi-purpose applications**

The OPTIMAX is an innovative table-top processor for film types up to a width of 35 cm (14"). The allrounder is used in various fields: human medicine including mammography, in veterinary medicine as well as for laboratory and graphic arts applications. The OPTIMAX is easy to service and presents itself with an attractive design.

#### **Uncompromising quality**

The one-piece casing and the exclusive use of high quality components guarantee the highest degree of reliability. It goes without saying that we supply the OPTIMAX with everything that is needed for operation like light protection cover, replenishment containers and connecting hoses.

#### **Best economic efficiency**

The integrated stand-by mode reduces your water and energy consumption. The film detection with a micro sensor enables intelligent and economic replenishment depending on the processed film surface. Because of that the OPTIMAX saves chemicals and provides a constant high quality of your image.

#### **No problem operation, perfect safety**

With the simple control panel the tanks can easily be filled during the warm-up cycle, and the developer temperature can be set. Acoustic and visual controls signal the readiness for processing of the next film. The anti-crystallisation and antioxidation programs in combination with the continuous circulation of the chemicals and the thermostat for the developer ensure constant high image quality. The processed film is available after 90 seconds. For special applications the processing time can be adjusted to 135 or 167 seconds (factory set. Please state preferred timing when placing order).

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**OPTIMAX. The multi-purpose precision processor with a wide range of accessories.**

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

*Extensive range of accessories: open base table, closed base cabinet and "through-the-wall-mounting"-kit.*

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*Coordinated product range: e.g. AIRCLEAN for fresh air in the darkroom or chemical mixer PROMIX A40.*

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## Technical Data OPTIMAX

Film format:	min. 10 x 10 cm (4" x 4"), max. 35 x ∞ cm (14" x ∞)	
Film type:	sheet and rollfilms	
Process time:	90/135/167 sec (factory set)	
Capacity at 90 sec process time:		
Film format (L x W in cm):	10 x 10 (4" x 4") (3 x parallel)	840 films/h
Film format (L x W in cm):	24 x 18 (9 1/2" x 7") (2 x parallel)	258 films/h
Film format (L x W in cm):	24 x 30 (9 1/2" x 12")	129 films/h
Film format (L x W in cm):	43 x 35 (17" x 14")	75 films/h
Developer temperature:	28 – 37°C (82 – 99°F)	
Dryer temperature:	65°C (149°F)	
Film detection:	microsensor	
Replenishment rate:	600 ml/m <sup>2</sup> (600 ml per 10 square feet)	
Replenishment containers:	optional 12 l or 25 l (3.2 or 6.6 gal)	
Anti-oxydation:	in 60 min intervals	
Anti-cristallisation:	in 20 min intervals	
Tank volumes D/F/W:	5/5/5 Litre (1.3/1.3/1.3 gal)	
Water consumption:	1.9 l/min (1/2 gal/min) during processing 0.0 l/min (0.0 gal/min) during stand-by	
Energy consumption:	1.4 kWh during processing 0.12 kWh during stand-by	
Dimensions (L/W/H):	77 x 59 x 42 cm (30" x 23" x 16 1/2")	
Weight (empty/filled):	35/50 kg (77/110 lb)	
Electrical connections:	220 – 240 V; 8.8 A; 50/60 Hz 110 – 120 V; 15 A (20 A plug); 50/60 Hz	
Water connection:	3/4"; 2 – 10 bar (1/4"; 29 – 145 psi)	

Ici est née une technique  
médicale innovante.  
*This is where innovative  
medical technology is produced.*