



Antepartum Fetal Monitor BT-300



Antepartum Fetal Monitor BT-300

- Compact and light
- Easy-to-use
- High sensitivity, probe-frequency 1 MHz
- Simultaneous external twin FHR with dual pulsed Doppler
- External uterine-contraction
- Fetal Movement
- Ergonomic designed event marker



Specifications

FHR

Dual Pulsed Doppler
 Ultrasound Frequency: 1 MHz
 Intensity: ≤ 10 mW/cm²
 FHR Range: 50 - 240 bpm
 FHR Accuracy: $\pm 2\%$ of Range

FM

Fetal Movement

UC

External Type
 Reference (zero) Control
 Measurement Range: 0-99 units

Printer

Thermal Array Type
 Print Speed: 1, 2, 3 cm/min
 Paper Feeding Function

Display

7-Segment LED
 3 Channels (FHR I, FHR II, UC)

Sound

Doppler Sound
 Alarm Sound
 Information Sound

Set-up

Alarm Upper / Lower Limit Value
 Print Speed, Print Contrast

Function

Mark Function
 FHR II Offset Function
 Auto Print Function

Power

Input: AC 100-240 V (50/60 Hz)
 Output: DC 16 V, 2.8 A
 Consumption: 80 VA

Dimension

191 (H) x 187 (L) x 201 (D) mm

Standard Accessory

Ultrasound Doppler Probes (2)
 UC Probe
 Event Marker Jack
 Print Paper
 Power Adaptor & Cord
 Ultrasound Gel
 Probe Belt
 Operation Manual

Warranty 2 Years