

# SonoTrax

## Handheld Fetal Heart Rate Doppler



SonoTrax Series Fetal Doppler can meet the routine examination requirements of obstetricians, midwives and pregnant women themselves.

SonoTrax Pro is the flagship of model with FHR digital LCD display and audio playback. SonoTrax Basic is a high performance model with FHR digital LCD display and SonoTrax Lite is for simple auscultation without LCD display.



SonoTrax Pro & Basic



SonoTrax Lite

### Features

- Ergonomic design, compact and light
- Easy-to-use, large LCD display and clean sound
- High sensitivity Interchangeable probe
- Accurate FHR detection
- Auto and Manual FHR Counting
- Optional rechargeable battery and charger
- Long time to continuous use (4 hours)
- Automatic power-off



Model	Interchangeable Probes	FHR Display	Built-in Audio Recorder	Backlight	Waterproof probe	Auto probe detection display	Battery status display	Low battery indication
SonoTrax Lite	Yes	No	No	No	Yes	No	No	Yes
SonoTrax Basic	Yes	Yes	No	Yes	Yes	Yes	Yes	Yes
SonoTrax Pro	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Accessory	Rechargeable Ni-H battery Re-charger Earphone Canvas Bag 9V Alkaline Battery							

